


Anti-Tac3 Rabbit Polyclonal Antibody For Wb, If-P, Elisa - Q9Uhf0

Item Number: CM0025737



Introduction

Rabbit polyclonal antibody targeting human TAC3 (Tachykinin-3), validated for WB, IF-P, ELISA. Cross-reacts with mouse and rat. Ideal for neuropeptide signaling and gonadal function research.

[Learn More](#)

Parameter	Value
Product Name	TAC3 Rabbit pAb
Remarks/Alias	NK3; NKB; HH10; NKNB; PRO1155; ZNEUROK1; LncZBTB39; TAC3
Species	Human
GeneID (Human)	6866
GeneID	6866
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UHF0
Protein Weight	13kDa
Shipping	Ice bag