

Anti-Tacc3 Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Y6A5

Item Number: CM0003750



Introduction

High-quality rabbit monoclonal TACC3 antibody for detection of endogenous TACC3 in human, mouse, and rat samples. Validated for Western blot, immunohistochemistry, immunofluorescence, and ELISA applications. Ideal for studying centrosome dynamics, mitosis, and cancer progression.

[Learn More](#)

Parameter	Value
Product Name	TACC3 Rabbit mAb
Remarks/Alias	ERIC1; Tacc4; ERIC-1; maskin; TACC3
Species	Human
Gene ID (Human)	10460
Gene ID	10460
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y6A5
Protein Weight	90kDa
Shipping	Ice bag