

Anti-Taco1 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Bsh4

Item Number: CM0017273

Introduction

Rabbit polyclonal antibody targeting human TACO1, a mitochondrial translational activator of cytochrome c oxidase I. Cross-reacts with mouse and rat. For WB, IHC-P, ELISA. Ideal for mitochondrial disease and respiratory chain research.

[Learn More](#)

Parameter	Value
Product Name	TACO1 Rabbit pAb
Remarks/Alias	CCDC44; MC4DN8; TACO1
Species	Human
GeneID (Human)	51204
GeneID	51204
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-297 of human TACO1 (NP_057444.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BSH4
Protein Weight	32kDa
Shipping	Ice bag