


Anti-Tacr3 Rabbit Polyclonal Antibody For Wb And Elisa - P29371

Item Number: CM0027970



Introduction

Rabbit polyclonal antibody against human TACR3 (Neuromedin-K receptor). Validated for Western blot and ELISA. Cross-reacts with mouse and rat. Recombinant immunogen (aa 356-465). Ideal for GPCR signaling research.

[Learn More](#)

Parameter	Value
Product Name	TACR3 Rabbit pAb
Remarks/Alias	NK3, NKR, HH11, NK3R, NK-3R, TAC3R, TAC3RL, TACR3
Species	Human
GeneID (Human)	6870
GeneID	6870
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 356-465 of human TACR3 (NP_001050.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P29371
Protein Weight	52kDa
Shipping	Ice bag