

Anti-Tada2B Polyclonal Antibody For Wb, Elisa - Q86Tj2

Item Number: CM0013314

Introduction

Rabbit polyclonal antibody targeting transcriptional adapter 2-beta (TADA2B), validated for WB and ELISA with cross-reactivity to human and mouse. Suitable for PAX5-dependent transcription and nuclear regulation studies.

[Learn More](#)

Parameter	Value
Product Name	TADA2B Rabbit pAb
Remarks/Alias	ADA2B; ADA2(beta); TADA2B
Species	Human
GeneID (human)	93624
GeneID	93624
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 150-322 of human TADA2B (NP_689506.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q86TJ2
Protein Weight	48kDa
Shipping	Ice bag