

Anti-Tada3 Polyclonal Antibody For Wb, If/Icc, Elisa - 075528

Item Number: CM0024880

Introduction

TADA3 rabbit polyclonal antibody validated for WB, IF/ICC, ELISA. Recognizes human, mouse, rat TADA3. Ideal for studying transcriptional regulation and histone acetyltransferase complexes.

[Learn More](#)

Parameter	Value
Product Name	TADA3 Rabbit pAb
Remarks/Alias	ADA3; NGG1; hADA3; STAF54; TADA3L; TADA3
Species	Human
GeneID (Human)	10474
GeneID	10474
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 270-369 of human TADA3 (NP_597814.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O75528
Protein Weight	49kDa
Shipping	Ice bag