

Anti-Taf1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P21675

Item Number: CM0014647



Introduction

Rabbit polyclonal antibody targeting human TAF1, the largest TFIID subunit. Validated for WB, IHC-P, ELISA; cross-reacts with mouse and rat. Ideal for transcriptional regulation research.

[Learn More](#)

Parameter	Value
Product Name	TAF1 Rabbit pAb
Remarks/Alias	OF; XDP; BA2R; CCG1; CCGS; DYT3; KAT4; P250; NSCL2; TAF2A; MRXS33; N-TAF1; TAFII250; DYT3/TAF1; TAFII-250; TAF(II)250; TAF1
Species	Human
GeneID (Human)	6872
GeneID	6872
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 990-1250 of human TAF1 (NP_620278.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P21675
Protein Weight	213kDa
Shipping	Ice bag