


Anti-Taf12 Polyclonal Antibody For Wb, Elisa - Q16514

Item Number: CM0024342



Introduction

Human TAF12 rabbit polyclonal antibody validated for Western blot and ELISA. Targets transcription initiation factor TFIID subunit 12, essential for RNA Pol II transcription. Ideal for transcription regulation studies.

[Learn More](#)

Parameter	Value
Product Name	TAF12 Rabbit pAb
Remarks/Alias	TAF2J; TAFII20; TAF12
Species	Human
GeneID (human)	6883
GeneID	6883
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-161 of human TAF12 (NP_001128690.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q16514
Protein Weight	18kDa
Shipping	Ice bag