

Anti-Taf15 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - Q92804

Item Number: CM0030436



Introduction

CamelBio supplies a high-quality TAF15 rabbit polyclonal antibody validated for WB, IHC-P, IF/ICC, IP, and ELISA. Essential for studying transcription initiation and RNA-binding mechanisms in human tissues.

[Learn More](#)

Parameter	Value
Product Name	TAF15 Rabbit pAb
Remarks/Alias	Npl3; RBP56; TAF2N; TAFII68; TAF15
Species	Human
GeneID (Human)	8148
GeneID	8148
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human TAF15 (NP_631961.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human
SWISS	Q92804
Protein Weight	62kDa
Shipping	Ice bag