

Anti-Taf1C Polyclonal Antibody For Wb, Elisa - Q15572

Item Number: CM0030055

Introduction

Rabbit polyclonal antibody targeting human TAF1C, suitable for Western blot and ELISA applications. Detects TAF1C from human, mouse, and rat samples. Ideal for studying RNA polymerase I transcription and nucleolar function.

[Learn More](#)

Parameter	Value
Product Name	TAF1C Rabbit pAb
Remarks/Alias	SL1; TAFI95; TAFI110; MGC:39976; TAF1C
Species	Human
GeneID (human)	9013
GeneID	9013
Immunogen	Recombinant protein, Recombinant fusion protein containing a sequence corresponding to amino acids 70-320 of human TAF1C (NP_001230085.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15572
Protein Weight	95kDa
Shipping	Ice bag