

Anti-Taf5 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q15542

Item Number: CM0029296



Introduction

Rabbit polyclonal antibody targeting human TAF5, a core TFIID subunit essential for RNA Polymerase II transcription. Validated for WB, IF/ICC, ELISA; cross-reacts with mouse and rat. Ideal for gene regulation and transcription studies.

[Learn More](#)

Parameter	Value
Product Name	TAF5 Rabbit pAb
Remarks/Alias	TAF2D; TAFII100; TAFII-100; TAF(II)100; TAF5
Species	Human
GeneID (Human)	6877
GeneID	6877
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 656-800 of human TAF5 (NP_008882.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15542
Protein Weight	87kDa
Shipping	Ice bag