

## Anti-Tafazzin / Taz Rabbit Polyclonal Antibody For Wb, Elisa - Q16635

Item Number: CM0016116

### Introduction

Tafazzin/TAZ rabbit polyclonal antibody for WB and ELISA. Detects human, mouse, rat endogenous Tafazzin. A mitochondrial acyltransferase essential for cardiolipin remodeling; suitable for Barth syndrome and cardiac muscle research.

[Learn More](#)

| Parameter        | Value                                                                                                                       |
|------------------|-----------------------------------------------------------------------------------------------------------------------------|
| Product Name     | Tafazzin / TAZ Rabbit pAb                                                                                                   |
| Remarks/Alias    | EFE; TAZ; BTHS; EFE2; G4.5; Taz1; CMD3A; LVNCX; Tafazzin / TAZ                                                              |
| Species          | Human                                                                                                                       |
| Gene ID (Human)  | 6901                                                                                                                        |
| Gene ID          | 6901                                                                                                                        |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 174-248 of human Tafazzin / TAZ (NP_851830.1) |
| Source           | Rabbit                                                                                                                      |
| Category         | Polyclonal Antibodies                                                                                                       |
| Application      | WB, ELISA                                                                                                                   |
| Cross Reactivity | Human, Mouse, Rat                                                                                                           |
| SWISS            | Q16635                                                                                                                      |
| Protein Weight   | 30kDa                                                                                                                       |
| Shipping         | Ice bag                                                                                                                     |