


Anti-Tak1 (Map3K7) Monoclonal Antibody For Wb, Ip, Elisa - O43318

Item Number: CM0003648



Introduction

Rabbit monoclonal antibody targeting human TAK1 (MAP3K7), validated for WB, IP, ELISA. Reacts with human and mouse. Useful for studying MAP kinase signaling, NF-kB activation, and apoptosis. Immunogen: synthetic peptide within amino acids 507-606 of human TAK1.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] TAK1 Rabbit mAb
Remarks/Alias	CSCF; FMD2; TAK1; MEKK7; TGF1a
Species	Human
GeneID (Human)	6885
GeneID	6885
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 507-606 of human TAK1 (O43318)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	O43318
Protein Weight	67 kDa
Shipping	Ice bag