


Anti-Tak1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O43318

Item Number: CM0017161



Introduction

Rabbit polyclonal antibody targeting human TAK1 (MAP3K7), validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat. Suitable for signal transduction, inflammation, and stress response research.

[Learn More](#)

Parameter	Value
Product Name	TAK1 Rabbit pAb
Remarks/Alias	CSCF; FMD2; TAK1; MEKK7; TGF1a
Species	Human
GeneID (Human)	6885
GeneID	6885
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 507-606 of human TAK1 (O43318)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43318
Protein Weight	67kDa
Shipping	Ice bag