

Anti-Tal2 Rabbit Polyclonal Antibody For Wb, Elisa - Q16559

Item Number: CM0013440

Introduction

Rabbit polyclonal antibody against TAL2 (bHLHa19), a basic helix-loop-helix transcription factor and proto-oncogene. Validated for WB and ELISA with cross-reactivity in rat. Supplied for leukemia and transcription research.

[Learn More](#)

Parameter	Value
Product Name	TAL2 Rabbit pAb
Remarks/Alias	TAL2
Species	Human
GeneID (human)	6887
GeneID	6887
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 20 to the C-terminus of human TAL2 (NP_005412.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	Q16559
Protein Weight	12kDa
Shipping	Ice bag