

Anti-Taldo1 Polyclonal Antibody For Wb, Elisa - P37837

Item Number: CM0013220



Introduction

Rabbit polyclonal antibody against human transaldolase (TALDO1) for WB and ELISA applications. Cross-reacts with mouse and rat. Supports pentose phosphate pathway and metabolic disorder research.

[Learn More](#)

Parameter	Value
Product Name	TALDO1 Rabbit pAb
Remarks/Alias	TAL; TALH; TAL-H; TALDOR; TALDO1
Species	Human
Gene ID (Human)	6888
Gene ID	6888
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-337 of human TALDO1 (NP_006746.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P37837
Protein Weight	38kDa
Shipping	Ice bag