

Anti-Taldo1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P37837

Item Number: CM0030936



Introduction

High-quality rabbit polyclonal antibody against human TALDO1, validated for WB, IF/ICC, and ELISA. Reacts with human and mouse TALDO1 (38 kDa), a key enzyme in the pentose phosphate pathway.

[Learn More](#)

Parameter	Value
Product Name	TALDO1 Rabbit pAb
Remarks/Alias	TAL; TALH; TAL-H; TALDOR; TALDO1
Species	Human
GeneID (human)	6888
GeneID	6888
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P37837
Protein Weight	38kDa
Shipping	Ice bag