

Anti-Talin-1 Polyclonal Antibody For Wb, If-P, Elisa - Q9Y490

Item Number: CM0031006



Introduction

Rabbit polyclonal antibody raised against human Talin-1. Validated for WB, IF-P, and ELISA, with cross-reactivity to human, mouse, and rat. Ideal for cytoskeleton and cell adhesion research.

[Learn More](#)

Parameter	Value
Product Name	TLN1 Rabbit pAb
Remarks/Alias	TLN; ILWEQ; talin-1; TLN1
Species	Human
Gene ID (Human)	7094
Gene ID	7094
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human TLN1 (NP_006280.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y490
Protein Weight	270kDa
Shipping	Ice bag