


Anti-Talin 2 Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - Q9Y4G6

Item Number: CM0003292



Introduction

Rabbit monoclonal Talin 2 antibody validated for Western blot, immunofluorescence, immunohistochemistry, and ELISA. Detects human, mouse, and rat Talin-2 (272 kDa). Ideal for cell adhesion and focal adhesion research.

[Learn More](#)

Parameter	Value
Product Name	Talin 2 Rabbit mAb
Remarks/Alias	ILWEQ, TLN2, Talin-2
Species	Human
Gene ID (Human)	83660
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 350-450 of human Talin 2 (Q9Y4G6)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y4G6
Protein Weight	272kDa
Shipping	Ice bag