

Anti-Tamm Horsfall Glycoprotein Polyclonal Antibody For Wb, If-P, Elisa - P07911

Item Number: CM0023982



Introduction

Rabbit polyclonal antibody against human Tamm Horsfall Glycoprotein (Uromodulin), validated for WB, IF-P, and ELISA. Cross-reacts with mouse and rat. Ideal for renal and urinary tract research.

[Learn More](#)

Parameter	Value
Product Name	Tamm Horsfall Glycoprotein Rabbit pAb
Remarks/Alias	THP; FJHN; HNFJ; THGP; HNFJ1; MCKD2; ADTKD1; ADMCKD2; Tamm Horsfall Glycoprotein
Species	Human
GeneID (Human)	7369
GeneID	7369
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 175-324 of human Tamm Horsfall Glycoprotein. (NP_003352.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P07911
Protein Weight	70kDa
Shipping	Ice bag