

Anti-Tamm41 Rabbit Polyclonal Antibody For Wb, Elisa - Q96Bw9

Item Number: CM0030267

Introduction

Rabbit polyclonal antibody against TAMM41 (phosphatidate cytidyltransferase, mitochondrial) for Western blot and ELISA. Reacts with human and mouse. Suitable for cardiolipin biosynthesis research.

[Learn More](#)

Parameter	Value
Product Name	TAMM41 Rabbit pAb
Remarks/Alias	RAM41, TAM41, C3orf31, COXPD56
Species	Human
GeneID (Human)	132001
GeneID	132001
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-316 of human TAMM41 (NP_620162.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q96BW9
Protein Weight	51kDa
Shipping	Ice bag