


Anti-Tank Rabbit Polyclonal Antibody For Wb, Elisa - Q92844

Item Number: CM0029680



Introduction

High-quality rabbit polyclonal antibody targeting human TANK (ITRAF), validated for Western blot and ELISA applications. Cross-reacts with human, mouse, and rat. Essential for NF- κ B signaling and innate immunity research.

[Learn More](#)

Parameter	Description
Product Name	TANK Rabbit pAb
Remarks/Alias	ITRAF; TRAF2; I-TRAF; TANK
Species	Human
GeneID (human)	10010
GeneID	10010
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-119 of human TANK (NP_597841.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92844
Protein Weight	48kDa
Shipping	Ice bag