

Anti-Taok1 Rabbit Monoclonal Antibody For Wb, Elisa - Q7L7X3

Item Number: CM0007388

Introduction

Anti-TAOK1 rabbit monoclonal antibody validated for WB and ELISA; reacts with human, mouse, and rat TAOK1 (116 kDa), a key kinase in p38 MAPK signaling, DNA damage response, and cytoskeletal regulation.

[Learn More](#)

Parameter	Value
Product Name	TAOK1 Rabbit mAb
Remarks/Alias	DDIB; PSK2; TAO1; KFC-B; MARKK; PSK-2; hKFC-B; hTAOK1; MAP3K16; TAOK1
Species	Human
GeneID (Human)	57551
GeneID	57551
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 900-985 of human TAOK1 (NP_065842.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q7L7X3
Protein Weight	116kDa
Shipping	Ice bag