

Anti-Taok2 Polyclonal Antibody For Wb And Elisa - Q9UI54

Item Number: CM0021570

Introduction

Rabbit polyclonal antibody against human TAOK2, a serine/threonine kinase involved in stress-activated MAPK cascades. Suitable for Western blot and ELISA. Recombinant immunogen; protein weight 138 kDa.

[Learn More](#)

Parameter	Value
Product Name	TAOK2 Rabbit pAb
Remarks/Alias	PSK; PSK1; TAO1; TAO2; MAP3K17; Tao2beta; PSK1-BETA; TAOK2
Species	Human
Gene ID (Human)	9344
Gene ID	9344
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 610-730 of human TAOK2 (NP_057235.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9UL54
Protein Weight	138kDa
Shipping	Ice bag