

Anti-Tap1 Monoclonal Antibody For Wb, Elisa - Q03518

Item Number: CM0004048



Introduction

High-quality rabbit monoclonal antibody targeting human TAP1 (Q03518) for Western blotting and ELISA. Cross-reacts with mouse TAP1; ideal for antigen processing and immunology studies.

[Learn More](#)

Parameter	Value
Product Name	TAP1 Rabbit mAb
Remarks/Alias	APT1; PSF1; ABC17; ABCB2; PSF-1; RING4; TAP1N; D6S114E; TAP1*0102N; TAP1
Species	Human
GeneID	6890
Immunogen	Recombinant protein corresponding to amino acids 449-748 of human TAP1 (NP_000584.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q03518
Protein Weight	81 kDa
Shipping	Ice bag