

Anti-Tap1 Polyclonal Antibody For Wb And Elisa - Q03518

Item Number: CM0030405



Introduction

Rabbit polyclonal antibody against human TAP1, validated for Western blot and ELISA. Cross-reacts with human and mouse. Ideal for studying MHC class I antigen processing and ABC transporter functions.

[Learn More](#)

Parameter	Value
Product Name	TAP1 Rabbit pAb
Remarks/Alias	APT1; PSF1; ABC17; ABCB2; PSF-1; RING4; TAP1N; D6S114E; TAP1*0102N; TAP1
Species	Human
Gene ID (Human)	6890
Gene ID	6890
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 524-808 of human TAP1 (NP_000584.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q03518
Protein Weight	81 kDa
Shipping	Ice bag