

## Anti-Tap2 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q03519

Item Number: CM0023823



### Introduction

TAP2 Rabbit Polyclonal Antibody validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Target is antigen peptide transporter 2, an ER membrane protein involved in MHC class I antigen presentation. Suitable for immunology studies.

[Learn More](#)

| Parameter        | Value                                                                                                                |
|------------------|----------------------------------------------------------------------------------------------------------------------|
| Product Name     | TAP2 Rabbit pAb                                                                                                      |
| Remarks/Alias    | APT2; PSF2; ABC18; ABCB3; PSF-2; RING11; D6S217E; TAP2                                                               |
| Species          | Human                                                                                                                |
| Gene ID (Human)  | 6891                                                                                                                 |
| Gene ID          | 6891                                                                                                                 |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 450-686 of human TAP2 (NP_001276972.1) |
| Source           | Rabbit                                                                                                               |
| Category         | Polyclonal Antibodies                                                                                                |
| Application      | WB, IHC-P, IF/ICC, ELISA                                                                                             |
| Cross Reactivity | Human, Mouse, Rat                                                                                                    |
| SWISS            | Q03519                                                                                                               |
| Protein Weight   | 76kDa                                                                                                                |
| Shipping         | Ice bag                                                                                                              |