

Anti-Tap2 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - Q03519

Item Number: CM0007292



Introduction

Rabbit monoclonal antibody targeting human TAP2, validated for WB, IF/ICC, and ELISA. Cross-reacts with mouse. Recognizes endogenous TAP2 at 76 kDa. Ideal for antigen processing and MHC class I presentation studies.

[Learn More](#)

Parameter	Value
Product Name	TAP2 Rabbit mAb
Remarks/Alias	APT2; PSF2; ABC18; ABCB3; PSF-2; RING11; D6S217E; TAP2
Species	Human
GeneID (human)	6891
GeneID	6891
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 450-686 of human TAP2 (NP_001276972.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q03519
Protein Weight	76kDa
Shipping	Ice bag