

Anti-Tapbp Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O15533

Item Number: CM0021354



Introduction

Rabbit polyclonal antibody targeting human tapasin (TAPBP) for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat; ideal for studying MHC class I assembly and peptide loading.

[Learn More](#)

Parameter	Value
Product Name	TAPBP Rabbit pAb
Remarks/Alias	TPN; TAPA; TPSN; NGS17; TAPBP
Species	Human
Gene ID (Human)	6892
Gene ID	6892
Immunogen	Recombinant protein containing a sequence corresponding to amino acids 30-260 of human TAPBP (NP_757345.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O15533
Protein Weight	48kDa
Shipping	Ice bag