

Anti-Tapp1 Rabbit Monoclonal Antibody For Wb, Elisa - Q9Hb21

Item Number: CM0002981



Introduction

High-specificity recombinant rabbit monoclonal antibody targeting human TAPP1 (PLEKHA1/TAPP-1). Validated for WB and ELISA, with cross-reactivity to mouse and rat. Ideal for phosphoinositide signaling and membrane trafficking research.

[Learn More](#)

Parameter	Value
Product Name	TAPP1 Rabbit mAb
Remarks/Alias	TAPP1
Species	Human
Gene ID (Human)	59338
Gene ID	59338
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 300-404 of human TAPP1 (Q9HB21).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9HB21
Protein Weight	46kDa
Shipping	Ice bag