

Anti-Tarbp1 Rabbit Polyclonal Antibody For Wb, Elisa - Q13395

Item Number: CM0020720



Introduction

High-quality rabbit polyclonal antibody targeting human TARBP1 (TRP-185), validated for WB and ELISA. Essential for studying tRNA methylation, HIV-1 TAR RNA binding, and related pathways.

[Learn More](#)

Parameter	Value
Product Name	TARBP1 Rabbit pAb
Remarks/Alias	TRM3; TRMT3; TRP185; TRP-185; TARBP1
Species	Human
GeneID (human)	6894
GeneID	6894
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1322-1621 of human TARBP1 (NP_005637.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q13395
Protein Weight	182kDa
Shipping	Ice bag