

Anti-Tarbp2 Rabbit Monoclonal Antibody For Wb, If, Ihc, Elisa - Q15633

Item Number: CM0011642



Introduction

Rabbit monoclonal TARBP2 antibody validated for WB, IF/ICC, IHC-P, ELISA. Reacts with human, mouse, rat. Essential for miRNA processing, HIV replication, and innate immunity research.

[Learn More](#)

Parameter	Value
Product Name	TARBP2 Rabbit mAb
Remarks/Alias	LOQS; TRBP; TRBP1; TRBP2; TARBP2
Species	Human
Gene ID (Human)	6895
Gene ID	6895
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 200-300 of human TARBP2 (Q15633).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15633
Protein Weight	39kDa
Shipping	Ice bag