

Anti-Tars Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P26639

Item Number: CM0026689



Introduction

Rabbit polyclonal antibody targeting human TARS (ThrRS). Validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Suitable for protein translation and aminoacyl-tRNA synthetase research.

[Learn More](#)

Parameter	Value
Product Name	TARS Rabbit pAb
Remarks/Alias	TARS; TTD7; ThrRS
Species	Human
Gene ID (Human)	6897
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-230 of human TARS (NP_689508.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P26639
Protein Weight	83kDa
Shipping	Ice bag