

Anti-Tars2 Polyclonal Antibody For Wb, Elisa - Q9Bw92

Item Number: CM0017014

Introduction

Rabbit polyclonal antibody targeting human TARS2, the mitochondrial threonine-tRNA ligase. Validated for Western blot and ELISA. Cross-reacts with human, mouse, and rat. Ideal for mitochondrial translation and disease studies.

[Learn More](#)

Parameter	Value
Product Name	TARS2 Rabbit pAb
Remarks/Alias	thrRS; TARSL1; COXPD21; TARS2
Species	Human
Gene ID	80222
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 100-400 of human TARS2 (NP_079426.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BW92
Protein Weight	81kDa
Shipping	Ice bag