

Anti-Tarsl2 Rabbit Polyclonal Antibody For Wb And Elisa - A2Rtx5

Item Number: CM0020135

Introduction

Anti-TARSL2 rabbit polyclonal antibody validated for WB and ELISA applications. Detects human, mouse, rat TARSL2 (ThrRS-L), a cytoplasmic threonine--tRNA ligase involved in protein biosynthesis. Suitable for IVD research.

[Learn More](#)

Parameter	Value
Product Name	TARSL2 Rabbit pAb
Remarks/Alias	TARSL2; ThrRS-L
Species	Human
GenelD (Human)	123283
GenelD	123283
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-170 of human TARSL2 (NP_689547.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	A2RTX5
Protein Weight	93kDa
Shipping	Ice bag