

Anti-Tas1R1 Rabbit Polyclonal Antibody For Wb And Elisa - Q7Rtx1

Item Number: CM0014599

Introduction

Rabbit polyclonal antibody targeting human TAS1R1 (GPR70), suitable for Western blot and ELISA. Detects umami taste receptor across human, mouse, and rat. Ideal for sensory transduction and GPCR studies.

[Learn More](#)

Parameter	Value
Product Name	TAS1R1 Rabbit pAb
Remarks/Alias	TR1, T1R1, GM148, GPR70, TAS1R1
Species	Human
GeneID (Human)	80835
GeneID	80835
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 330-567 of human TAS1R1 (NP_619642.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q7RTX1
Protein Weight	93kDa
Shipping	Ice bag