

Anti-Tas1R2 Rabbit Polyclonal Antibody For Wb And Elisa - Q8Te23

Item Number: CM0023881



Introduction

Rabbit polyclonal antibody against human TAS1R2 (sweet taste receptor T1R2), validated for Western blot and ELISA. Supports taste and sensory transduction studies.

[Learn More](#)

Parameter	Value
Product Name	TAS1R2 Rabbit pAb
Remarks/Alias	TR2; T1R2; GPR71; TAS1R2
Species	Human
Gene ID (Human)	80834
Gene ID	80834
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q8TE23
Protein Weight	95kDa
Shipping	Ice bag