


Anti-Tas1R3 Polyclonal Antibody For Wb, Elisa - Q7Rtx0

Item Number: CM0016088



Introduction

Rabbit polyclonal antibody against human TAS1R3, suitable for WB and ELISA. Cross-reacts with mouse and rat. Ideal for taste receptor research, umami and sweet taste signaling studies.

[Learn More](#)

Parameter	Value
Product Name	TAS1R3 Rabbit pAb
Remarks/Alias	T1R3; TAS1R3
Species	Human
Gene ID (Human)	83756
Gene ID	83756
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 400-570 of human TAS1R3 (NP_689414.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q7RTX0
Protein Weight	93kDa
Shipping	Ice bag