

Anti-Tas2R10 Polyclonal Antibody For Wb, Elisa - Q9Nyw0

Item Number: CM0017847

Introduction

Rabbit polyclonal antibody against human taste receptor TAS2R10 (TRB2/T2R10). Validated for WB and ELISA, with cross-reactivity to mouse and rat. For bitter taste perception and gustducin signaling research.

[Learn More](#)

Parameter	Value
Product Name	TAS2R10 Rabbit pAb
Remarks/Alias	TRB2; T2R10; TAS2R10
Species	Human
Gene ID (Human)	50839
Gene ID	50839
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 100-200 of human TAS2R10 (NP_076410.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9Nyw0
Protein Weight	35kDa
Shipping	Ice bag