

Anti-Tas2R38 Rabbit Polyclonal Antibody For Wb, Elisa - P59533

Item Number: CM0016517

Introduction

Rabbit polyclonal antibody against human TAS2R38, a bitter taste receptor. Suitable for WB and ELISA, with cross-reactivity to rat. Ideal for investigating bitter taste perception and GPCR signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	TAS2R38 Rabbit pAb
Remarks/Alias	PTC, T2R38, T2R61, THIOT, TAS2R38
Species	Human
GeneID (Human)	5726
GeneID	5726
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 150-250 of human TAS2R38 (NP_789787.4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	P59533
Protein Weight	38kDa
Shipping	Ice bag