


Anti-Tasl Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - Q9Hai6

Item Number: CM0012381



Introduction

Anti-TASL rabbit monoclonal antibody validated for Western blot, IF/ICC, and ELISA. Targets human TASL (CXorf21; 34 kDa), a lysosomal adapter in TLR7/8/9-IRF5 innate immune signaling. Suitable for immunology research.

[Learn More](#)

Parameter	Specification
Product Name	TASL Rabbit mAb
Remarks/Alias	CXorf21
Species	Human
GeneID (Human)	80231
GeneID	80231
Immunogen	Synthetic peptide corresponding to amino acids 211-301 of human TASL (NP_079435.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q9HAI6
Protein Weight	34 kDa
Shipping	Ice bag