


Anti-Tasl Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Hai6

Item Number: CM0024695



Introduction

Rabbit polyclonal antibody against human TASL (CXorf21), validated for Western blot, IHC-P, and ELISA. Cross-reacts with mouse. Suitable for innate immunity and TLR signaling research.

[Learn More](#)

Parameter	Value
Product Name	TASL Rabbit pAb
Remarks/Alias	CXorf21; TASL
Species	Human
GeneID (Human)	80231
GeneID	80231
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 96-196 of human TASL (NP_079435.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9HAI6
Protein Weight	34kDa
Shipping	Ice bag