


Anti-Tat Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P17735

Item Number: CM0025313



Introduction

Rabbit polyclonal antibody against human TAT (tyrosine aminotransferase), validated for WB, IHC-P, ELISA. Cross-reacts with human, mouse, rat. 50kDa. Ideal for tyrosine metabolism studies.

[Learn More](#)

Parameter	Value
Product Name	TAT Rabbit pAb
Remarks/Alias	TAT
Species	Human
GeneID (human)	6898
GeneID	6898
Immunogen	Recombinant protein / Recombinant fusion protein containing a sequence corresponding to amino acids 1-270 of human TAT (NP_000344.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P17735
Protein Weight	50kDa
Shipping	Ice bag