


Anti-Tatdn1 Polyclonal Antibody For Wb And Elisa - Q6P1N9

Item Number: CM0016942



Introduction

Rabbit polyclonal antibody against human TATDN1, a nuclear deoxyribonuclease with exonuclease/endonuclease activities. Validated for WB, ELISA. Useful for DNA repair and chromosomal segregation research.

[Learn More](#)

Parameter	Value
Product Name	TATDN1 Rabbit pAb
Remarks/Alias	CDA11; TATDN1
Species	Human
GeneID (Human)	83940
GeneID	83940
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-297 of human TATDN1 (NP_114415.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q6P1N9
Protein Weight	34kDa
Shipping	Ice bag