

Anti-Tau Monoclonal Antibody For Wb - P10636

Item Number: CM0006279



Introduction

Rabbit monoclonal antibody targeting human Tau (MAPT, 79 kDa). Validated for WB, IP, IF, IHC-P, ELISA, with cross-reactivity to mouse and rat. Ideal for Alzheimer's and neurodegeneration research.

[Learn More](#)

Parameter	Value
Product Name	Tau Rabbit mAb
Remarks/Alias	TAU; MSTD; PPND; DDPAC; MAPTL; MTBT1; MTBT2; tau-40; FTDP-17; PPP1R103; Tau-PHF6; Tau
Species	Human
GeneID	4137
Immunogen	Synthetic Peptide - A synthetic peptide corresponding to a sequence within amino acids 1-100 of human MAPT (NM_005910.6)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF-F, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P10636
Protein Weight	79kDa
Shipping	Ice bag