

Anti-Tau Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P10636

Item Number: CM0016298



Introduction

Rabbit polyclonal antibody against human Tau protein (aa 658-758). For WB, IHC-P, IF/ICC, ELISA; cross-reacts with human, mouse, rat. Tau is a microtubule-associated protein crucial for neuronal polarity and stability.

[Learn More](#)

Parameter	Value
Product Name	Tau Rabbit pAb
Remarks/Alias	TAU; MSTD; PPND; DDPAC; MAPTL; MTBT1; MTBT2; tau-40; FTDP-17; PPP1R103; Tau-PHF6; Tau
Species	Human
Gene ID (Human)	4137
Gene ID	4137
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 658-758 of human Tau (NP_058519.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P10636
Protein Weight	79kDa
Shipping	Ice bag