

## Anti-Tau Rabbit Polyclonal Antibody For Wb, If/Icc, If-P, Ihc-P, Elisa - P10636

Item Number: CM0017268



### Introduction

Rabbit polyclonal antibody targeting human Tau protein (MAPT). Validated for WB, IF/ICC, IF-P, IHC-P, ELISA with cross-reactivity to mouse and rat. Ideal for studying microtubule dynamics and neurodegeneration research.

[Learn More](#)

| Parameter        | Value                                                                                               |
|------------------|-----------------------------------------------------------------------------------------------------|
| Product Name     | Tau Rabbit pAb                                                                                      |
| Remarks/Alias    | TAU; MSTD; PPND; DDPAC; MAPTL; MTBT1; MTBT2; tau-40; FTDP-17; PPP1R103; Tau-PHF6; Tau               |
| Species          | Human                                                                                               |
| GeneID (Human)   | 4137                                                                                                |
| GeneID           | 4137                                                                                                |
| Immunogen        | Synthetic peptide corresponding to a sequence within amino acids 500-758 of human Tau (NP_058519.3) |
| Source           | Rabbit                                                                                              |
| Category         | Polyclonal Antibodies                                                                               |
| Application      | WB, IF/ICC, IF-P, IHC-P, ELISA                                                                      |
| Cross Reactivity | Human, Mouse, Rat                                                                                   |
| SWISS            | P10636                                                                                              |
| Protein Weight   | 79kDa                                                                                               |
| Shipping         | Ice bag                                                                                             |