


## Anti-Tax1BP3 Polyclonal Antibody For Wb, Elisa - O14907

Item Number: CM0026341



### Introduction

Anti-TAX1BP3 Polyclonal Antibody (Rabbit) for WB and ELISA applications, reactive with Human and Mouse. Targets Tax1-binding protein 3, a Wnt signaling inhibitor. Ideal for PDZ domain and Rho signaling research.

[Learn More](#)

| Parameter        | Value                                                                                                               |
|------------------|---------------------------------------------------------------------------------------------------------------------|
| Product Name     | TAX1BP3 Rabbit pAb                                                                                                  |
| Remarks/Alias    | TIP1; TIP-1; TAX1BP3                                                                                                |
| Species          | Human                                                                                                               |
| GeneID (human)   | 30851                                                                                                               |
| GeneID           | 30851                                                                                                               |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-124 of human TAX1BP3 (NP_055419.1). |
| Source           | Rabbit                                                                                                              |
| Category         | Polyclonal Antibodies                                                                                               |
| Application      | WB, ELISA                                                                                                           |
| Cross Reactivity | Human, Mouse                                                                                                        |
| SWISS            | O14907                                                                                                              |
| Protein Weight   | 14kDa                                                                                                               |
| Shipping         | Ice bag                                                                                                             |