


Anti-Tax1BP3 Rabbit Monoclonal Antibody For Wb And Elisa - O14907

Item Number: CM0003554



Introduction

Rabbit monoclonal antibody targeting TAX1BP3 (TIP-1), validated for WB and ELISA. Reacts with human, mouse, and rat TAX1BP3. Suitable for studies on Wnt signaling, PDZ interactions, and viral oncogenesis.

[Learn More](#)

Parameter	Value
Product Name	TAX1BP3 Rabbit mAb
Remarks/Alias	TIP1; TIP-1; TAX1BP3
Species	Human
GeneID (Human)	30851
GeneID	30851
Immunogen	Synthetic Peptide; A synthetic peptide corresponding to a sequence within amino acids 1-100 of human TAX1BP3 (NP_055419.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O14907
Protein Weight	14kDa
Shipping	Ice bag