


Anti-Taz Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Gzv5

Item Number: CM0027978



Introduction

Anti-TAZ rabbit polyclonal antibody validated for WB, IHC-P, IF/ICC, ELISA applications. Reacts with human, mouse, and rat samples. Ideal for Hippo signaling, transcription regulation, and cancer research studies.

[Learn More](#)

Parameter	Value
Product Name	TAZ Rabbit pAb
Remarks/Alias	TAZ
Species	Human
Gene ID	25937
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-230 of human TAZ (NP_056287.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9GZV5
Protein Weight	44kDa
Shipping	Ice bag